


Anti-Mark1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9POL2

Item Number: CM0027156



Introduction

MARK1 polyclonal antibody validated for WB, IHC-P, ELISA. Detects human and mouse MARK1, a serine/threonine kinase regulating cell polarity and Wnt signaling. Ideal for neuroscience and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MARK1 Rabbit pAb
Remarks/Alias	MARK; Par1c; Par-1c; MARK1
Species	Human
GeneID (Human)	4139
GeneID	4139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 560-700 of human MARK1 (NP_061120.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9POL2
Protein Weight	89kDa
Shipping	Ice bag