

## Anti-Mark2 Rabbit Monoclonal Antibody For Wb And Elisa - Q7Kzi7

Item Number: CM0008334



### Introduction

High specificity MARK2 (EMK1/PAR-1) rabbit monoclonal antibody validated for WB and ELISA, cross-reactive with human, mouse, and rat. Ideal for studying cell polarity, microtubule dynamics, and neurobiology.

[Learn More](#)

| Parameter        | Value                                                                                            |
|------------------|--------------------------------------------------------------------------------------------------|
| Product Name     | MARK2 Rabbit mAb                                                                                 |
| Remarks/Alias    | EMK1; EMK-1; PAR-1; Par1b; Par-1b; MARK2                                                         |
| Species          | Human                                                                                            |
| Gene ID          | 2011                                                                                             |
| Immunogen        | Synthetic Peptide corresponding to a sequence within amino acids 689-788 of human MARK2 (Q7KZI7) |
| Source           | Rabbit                                                                                           |
| Category         | Monoclonal Antibodies                                                                            |
| Application      | WB, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                |
| SWISS            | Q7KZI7                                                                                           |
| Protein Weight   | 88kDa                                                                                            |
| Shipping         | Ice bag                                                                                          |