


Anti-Mark2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7Kzi7

Item Number: CM0021910



Introduction

MARK2 Rabbit pAb validated for WB, IF/ICC, ELISA with human, mouse, rat cross-reactivity. Ideal for cell polarity, microtubule dynamics, and neuronal migration research.

[Learn More](#)

Parameter	Value
Product Name	MARK2 Rabbit pAb
Remarks/Alias	EMK1; EMK-1; PAR-1; Par1b; Par-1b; MARK2
Species	Human
GeneID (Human)	2011
GeneID	2011
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-455 of human MARK2 (Q7KZI7).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7KZI7
Protein Weight	88kDa
Shipping	Ice bag