

Anti-Mark3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P27448

Item Number: CM0019155



Introduction

Rabbit polyclonal antibody targeting human MARK3 kinase, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. MARK3 phosphorylates microtubule-associated proteins and regulates Hippo signaling and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	MARK3 Rabbit pAb
Remarks/Alias	KP78; VIPB; CTAK1; PAR1A; Par-1a; MARK3
Species	Human
Gene ID (Human)	4140
Gene ID	4140
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MARK3 (NP_001122390.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27448
Protein Weight	84kDa
Shipping	Ice bag