

Anti-Marveld1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bsk0

Item Number: CM0019342

Introduction

Rabbit polyclonal antibody against human MARVELD1 (Gene ID: 83742) for WB and ELISA applications. Validated in human, mouse, and rat. Useful for studying cell cycle regulation and tumor suppression mechanisms.

[Learn More](#)

Parameter	Value
Product Name	MARVELD1 Rabbit pAb
Remarks/Alias	GB14; MARVD1; MRVLDC1; bA548K23.8; MARVELD1
Species	Human
Gene ID (Human)	83742
Gene ID	83742
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-50 of human MARVELD1 (NP_113672.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSK0
Protein Weight	19kDa
Shipping	Ice bag