


Anti-Marveld2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N4S9

Item Number: CM0025962



Introduction

Rabbit polyclonal antibody targeting MARVELD2 (tricellulin), validated for WB and ELISA. Reacts with human, mouse, and rat. Essential for tight junction and hearing loss research.

[Learn More](#)

Parameter	Value
Product Name	MARVELD2 Rabbit pAb
Remarks/Alias	Tric; DFNB49; MARVD2; MRVLDC2; MARVELD2
Species	Human
GeneID	153562
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human MARVELD2 (NP_001033692.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4S9
Protein Weight	64kDa
Shipping	Ice bag