

Anti-Mas1L Rabbit Polyclonal Antibody For Wb, Elisa - P35410

Item Number: CM0013969

Introduction

Rabbit polyclonal antibody targeting human MAS1L (MRG, MAS1-like), validated for WB and ELISA. Ideal for membrane receptor studies and G protein-coupled receptor research.

[Learn More](#)

Parameter	Value
Product Name	MAS1L Rabbit pAb
Remarks/Alias	MRG; MAS-L; dJ994E9.2; MAS1L
Species	Human
GeneID (Human)	116511
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human MAS1L (NP_443199.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35410
Protein Weight	42kDa
Shipping	Ice bag