

Anti-Mast Cell Chymase (Cma1) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P23946

Item Number: CM0002855



Introduction

High-quality rabbit monoclonal antibody against human Mast Cell Chymase (CMA1). Ideal for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Validated for chymase protease studies in mast cell biology.

[Learn More](#)

Parameter	Value
Product Name	Mast Cell Chymase (CMA1) Rabbit mAb
Remarks/Alias	CYH; MCT1; chymase; Mast Cell Chymase (CMA1)
Species	Human
GeneID (human)	1215
GeneID	1215
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 148-247 of human Mast Cell Chymase (CMA1) (CMA1) (P23946).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23946
Protein Weight	27kDa
Shipping	Ice bag