


Anti-Mast Cell Tryptase (Tpsb2) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15661

Item Number: CM0002881



Introduction

Rabbit monoclonal antibody targeting Mast Cell Tryptase (TPSB2/TPSAB1). Validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Ideal for mast cell biology and innate immunity research.

[Learn More](#)

Parameter	Description
Product Name	Mast Cell Tryptase (TPSB2) Rabbit mAb
Remarks/Alias	TPS1; TPS2; TPSB1; TPSB2; Tryptase-2; Mast Cell Tryptase (TPSB2)
Species	Human
Gene ID (Human)	7177
Gene ID	64499
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-275 of human Mast Cell Tryptase (TPSB2) (P20231)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15661
Protein Weight	31kDa
Shipping	Ice bag