

Anti-Mat2A Polyclonal Antibody For Wb, Elisa - P31153

Item Number: CM0029885



Introduction

MAT2A rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets S-adenosylmethionine synthase isoform type-2, a key enzyme in one-carbon metabolism. Ideal for research in methylation pathways.

[Learn More](#)

Parameter	Value
Product Name	MAT2A Rabbit pAb
Remarks/Alias	MATA2; MATII; SAMS2; MAT2A
Species	Human
GeneID	4144
Immunogen	Recombinant fusion protein containing amino acids 1-395 of human MAT2A (NP_005902.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31153
Protein Weight	44kDa
Shipping	Ice bag