

Anti-Mat2A Rabbit Monoclonal Antibody For Wb, Elisa - P31153

Item Number: CM0003556



Introduction

Rabbit monoclonal antibody against human MAT2A, validated for WB and ELISA. Reacts with human, mouse, and rat. Targets S-adenosylmethionine synthase isoform type-2, key in one-carbon metabolism.

[Learn More](#)

Parameter	Value
Product Name	MAT2A Rabbit mAb
Remarks/Alias	MATA2; MATII; SAMS2; MAT2A
Species	Human
GeneID (Human)	4144
GeneID	4144
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAT2A (P31153).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31153
Protein Weight	44kDa
Shipping	Ice bag