


Anti-Matk Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P42679

Item Number: CM0029270



Introduction

High-quality anti-MATK rabbit polyclonal antibody for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat MATK. Ideal for hematopoietic signal transduction and kinase regulation studies.

[Learn More](#)

Parameter	Value
Product Name	MATK Rabbit pAb
Remarks/Alias	CHK; CTK; HYL; Lsk; HYLTK; HHYLTK; MATK
Species	Human
Gene ID (Human)	4145
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 352-466 of human MATK (NP_647611.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42679
Protein Weight	56kDa
Shipping	Ice bag