

Anti-Matn3 Rabbit Polyclonal Antibody For Wb, Elisa - 015232

Item Number: CM0027276

Introduction

High-quality rabbit polyclonal antibody targeting human MATN3, validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cartilage extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	MATN3 Rabbit pAb
Remarks/Alias	HOA; OS2; EDM5; DIPOA; OADIP; SEMDBCD; MATN3
Species	Human
Gene ID	4148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 277-486 of human MATN3 (NP_002372.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15232
Protein Weight	53kDa
Shipping	Ice bag