

Anti-Matn4 Rabbit Polyclonal Antibody For Wb And Elisa - O95460

Item Number: CM0027063

Introduction

Rabbit polyclonal antibody targeting human MATN4 (matrilin-4), a cartilage extracellular matrix protein. Validated for WB, ELISA; cross-reacts with human, mouse, rat. Suitable for ECM and skeletal development research.

[Learn More](#)

Parameter	Value
Product Name	MATN4 Rabbit pAb
Remarks/Alias	MATN4
Species	Human
Gene ID (Human)	8785
Gene ID	8785
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 221-540 of human MATN4 (NP_085080.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95460
Protein Weight	68kDa
Shipping	Ice bag