

Anti-Matrilin-2 (Matn2) Rabbit Polyclonal Antibody For Wb And Elisa - O00339

Item Number: CM0021006

Introduction

Rabbit polyclonal antibody targeting human Matrilin-2 (MATN2), validated for Western blot (WB) and ELISA. Detects a 107 kDa protein; cross-reacts with human and rat. Ideal for extracellular matrix assembly research.

[Learn More](#)

Parameter	Value
Product Name	MATN2 Rabbit pAb
Remarks/Alias	MATN2
Species	Human
GeneID (Human)	4147
GeneID	4147
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 800 to the C-terminus of human MATN2 (NP_001304677.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O00339
Protein Weight	107kDa
Shipping	Ice bag