

Anti-Matrin 3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43243

Item Number: CM0002302

Introduction

This Matrin 3 rabbit monoclonal antibody specifically detects human Matrin 3 (P43243) and cross-reacts with mouse and rat. Suitable for WB, IHC-P, IF/ICC, and ELISA, target size ~95 kDa.

[Learn More](#)

Parameter	Value
Product Name	Matrin 3 Rabbit mAb
Remarks/Alias	MPD2; ALS21; VCPDM; Matrin 3
Species	Human
Gene ID (Human)	9782
Gene ID	9782
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 700-800 of human Matrin 3 (P43243).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43243
Protein Weight	95kDa
Shipping	Ice bag