

Anti-Matrin 3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43243

Item Number: CM0028922

Introduction

Matrin 3 Rabbit Polyclonal Antibody (P43243) validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for studying nuclear matrix, transcription, RNA binding, and innate immunity signaling.

[Learn More](#)

Parameter	Value
Product Name	Matrin 3 Rabbit pAb
Remarks/Alias	MPD2; ALS21; VCPDM; Matrin 3
Species	Human
GenelD (Human)	9782
GenelD	9782
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 747-847 of human Matrin 3 (NP_001181884.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43243
Protein Weight	95kDa
Shipping	Ice bag