

Anti-Matrix Metalloproteinase-11 (Mmp11/Stromelysin-3) Rabbit Monoclonal Antibody For Wb, Elisa - P24347

Item Number: CM0002042

Introduction

Recombinant rabbit monoclonal antibody against human MMP11, validated for WB and ELISA. Detects MMP11 in human, mouse, and rat samples. Suitable for cancer and stromal biology research.

[Learn More](#)

Parameter	Value
Product Name	MMP11 Rabbit mAb
Remarks/Alias	ST3; SL-3; STMY3; MMP11
Species	Human
GeneID (Human)	4320
GeneID	4320
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 389-488 of human MMP11 (P24347).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24347
Protein Weight	55kDa
Shipping	Ice bag