


Anti-Mbd2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ubb5

Item Number: CM0008108



Introduction

Anti-MBD2 Rabbit Monoclonal Antibody detects human and mouse MBD2 by WB and ELISA. Ideal for research on DNA methylation, chromatin remodeling, and gene silencing in cancer and epigenetics.

[Learn More](#)

Parameter	Value
Product Name	MBD2 Rabbit mAb
Remarks/Alias	DMTase; NY-CO-41; MBD2
Species	Human
Gene ID (Human)	8932
Gene ID	8932
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human MBD2 (Q9UBB5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UBB5
Protein Weight	43kDa
Shipping	Ice bag