


Anti-Mbd3 Monoclonal Antibody For Wb, Ihc-P, Elisa - O95983

Item Number: CM0006208



Introduction

High-quality rabbit monoclonal antibody against human MBD3. Validated for Western blot, IHC-P, and ELISA. Recognizes human, mouse, and rat MBD3. Suitable for chromatin remodeling and gene silencing studies.

[Learn More](#)

Parameter	Value
Product Name	MBD3 Rabbit mAb
Remarks/Alias	MBD3; methyl-CpG-binding domain protein 3
Species	Human
GeneID (Human)	53615
GeneID	53615
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human MBD3 (O95983)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95983
Protein Weight	33kDa
Shipping	Ice bag