

Anti-Mbl2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P11226

Item Number: CM0025177



Introduction

Rabbit polyclonal antibody targeting human MBL2 (mannose-binding lectin), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for studying the lectin complement pathway in innate immunity research. - P11226

[Learn More](#)

Parameter	Value
Product Name	MBL2 Rabbit pAb
Remarks/Alias	MBL; MBP; MBP1; MBPD; MBL2D; MBP-C; COLEC1; HSMBPC; MBL2
Species	Human
Gene ID (Human)	4153
Gene ID	4153
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-248 of human MBL2 (NP_000233.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11226
Protein Weight	26kDa
Shipping	Ice bag