

Anti-Mbnl1 Polyclonal Antibody For Wb And Elisa - Q9Nr56

Item Number: CM0026895



Introduction

Rabbit polyclonal antibody against human MBNL1, validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying alternative splicing regulation in muscle and cardiac biology.

[Learn More](#)

Parameter	Value
Product Name	MBNL1 Rabbit pAb
Remarks/Alias	EXP; MBNL
Species	Human
Gene ID (Human)	4154
Gene ID	4154
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 280-380 of human MBNL1 (NP_997176.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NR56
Protein Weight	42kDa
Shipping	Ice bag