

## Anti-Mbnl1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Nr56

Item Number: CM0004273



### Introduction

High-quality MBNL1 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA applications with human specificity. Ideal for research on alternative splicing regulation mediated by muscleblind-like protein 1, including cardiac troponin-T and insulin receptor pre-mRNA processing.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | MBNL1 Rabbit mAb                                                                                    |
| Remarks/Alias    | EXP; MBNL; MBNL1                                                                                    |
| Species          | Human                                                                                               |
| Gene ID          | 4154                                                                                                |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 289-388 of human MBNL1 (Q9NR56). |
| Source           | Rabbit                                                                                              |
| Category         | Monoclonal Antibody                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                   |
| Cross Reactivity | Human                                                                                               |
| SWISS            | Q9NR56                                                                                              |
| Protein Weight   | 42kDa                                                                                               |
| Shipping         | Ice bag                                                                                             |