

Anti-Mbnl1 Rabbit Polyclonal Antibody For Wb - Q9Nr56

Item Number: CM0030804



Introduction

Rabbit polyclonal antibody targeting human MBNL1, suitable for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for alternative splicing and muscle biology studies.

[Learn More](#)

Parameter	Value
Product Name	MBNL1 Rabbit pAb
Remarks/Alias	EXP; MBNL; MBNL1
Species	Human
GeneID (Human)	4154
GeneID	4154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-382 of human MBNL1 (NP_997176.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR56
Protein Weight	42kDa
Shipping	Ice bag