

Anti-Mbtps1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14703

Item Number: CM0023050



Introduction

Rabbit polyclonal antibody against human MBTPS1 (S1P/SKI-1), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for research on cholesterol metabolism, ER stress, lysosomal biogenesis, and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	MBTPS1 Rabbit pAb
Remarks/Alias	S1P; PCSK8; SEDKF; SKI-1; MBTPS1
Species	Human
GeneID (Human)	8720
GeneID	8720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 803-1052 of human MBTPS1 (NP_003782.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14703
Protein Weight	118kDa
Shipping	Ice bag