

Anti-Mbtps1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14703

Item Number: CM0004363



Introduction

Rabbit monoclonal antibody against human MBTPS1 (S1P/PCSK8), validated in WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Suitable for studying SREBP cleavage, lysosomal biogenesis, and ER stress pathways.

[Learn More](#)

Parameter	Value
Product Name	MBTPS1 Rabbit PolymAb®
Remarks/Alias	S1P; PCSK8; SEDKF; SKI-1
Species	Human
Gene ID (Human)	8720
Gene ID	8720
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 953-1052 of human MBTPS1 (NP_003782.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14703
Protein Weight	118kDa
Shipping	Ice bag