

Anti-Mc1 Receptor Monoclonal Antibody For Wb, Elisa - Q01726

Item Number: CM0002968



Introduction

Rabbit monoclonal antibody against human MC1 Receptor (MSH-R) for WB and ELISA. Cross-reacts with mouse and rat. Validated for research on melanocortin signaling and pigmentation.

[Learn More](#)

Parameter	Value
Product Name	MC1 Receptor Rabbit mAb
Remarks/Alias	CMM5; MSH-R; SHEP2; MC1 Receptor
Species	Human
GeneID (Human)	4157
GeneID	4157
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human MC1 Receptor (Q01726)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01726
Protein Weight	35kDa
Shipping	Ice bag