

Anti-Mc1 Receptor Polyclonal Antibody For Wb, Elisa - Q01726

Item Number: CM0014274



Introduction

Anti-MC1 Receptor Rabbit pAb validated for WB and ELISA. Cross-reacts with human, mouse, rat. Targets melanocyte-stimulating hormone receptor involved in melanogenesis. Ideal for GPCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	MC1 Receptor Rabbit pAb
Remarks/Alias	CMM5; MSH-R; SHEP2; MC1 Receptor
Species	Human
Gene ID (Human)	4157
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-317 of human MC1 Receptor (NP_002377.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01726
Protein Weight	35kDa
Shipping	Ice bag