


Anti-Mc5 Receptor Rabbit Monoclonal Antibody For Wb, Elisa - P33032

Item Number: CM0011347



Introduction

Rabbit monoclonal antibody targeting human MC5 Receptor (MC5R, MC-2) for WB and ELISA. Cross-reacts with mouse and rat. Ideal for GPCR signaling and melanocortin research. Predicted MW 37 kDa.

[Learn More](#)

Parameter	Value
Product Name	MC5 Receptor Rabbit mAb
Remarks/Alias	MC2; MC5 Receptor
Species	Human
GeneID	4161
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MC5 Receptor (P33032).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33032
Protein Weight	37kDa
Shipping	Ice bag