


Anti-Mcat Rabbit Polyclonal Antibody For Wb And Elisa - Q8lvs2

Item Number: CM0017562



Introduction

MCAT Rabbit pAb for Western blot (WB) and ELISA applications; targets human malonyl-CoA-acyl carrier protein transacylase, mitochondrial; recombinant immunogen, 43kDa. Ideal for fatty acid biosynthesis and mitochondrial biology research.

[Learn More](#)

Parameter	Value
Product Name	MCAT Rabbit pAb
Remarks/Alias	MT; MCT; MCT1; fabD; NET62; FASN2C; MCAT
Species	Human
GeneID (human)	27349
GeneID	27349
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 191-390 of human MCAT (NP_775738.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IVS2
Protein Weight	43kDa
Shipping	Ice bag