


Anti-Mccc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Rq3

Item Number: CM0014812



Introduction

This rabbit polyclonal antibody targets human MCCC1, validated for WB, IHC-P, and ELISA applications. It cross-reacts with mouse and rat, supporting research in leucine metabolism and related disorders.

[Learn More](#)

Parameter	Value
Product Name	MCCC1 Rabbit pAb
Remarks/Alias	MCCA; MCC-B; MCCCalpha; MCCC1
Species	Human
Gene ID (Human)	56922
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-725 of human MCCC1 (NP_064551.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RQ3
Protein Weight	80kDa
Shipping	Ice bag