


Anti-Mccc2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hcc0

Item Number: CM0017390



Introduction

Rabbit polyclonal antibody against MCCC2 (MCCB), validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Suitable for studying methylcrotonoyl-CoA carboxylase in leucine catabolism.

[Learn More](#)

Parameter	Value
Product Name	MCCC2 Rabbit pAb
Remarks/Alias	MCCB; MCCCbeta; MCCC2
Species	Human
GeneID (Human)	64087
GeneID	64087
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 284-563 of human MCCC2 (NP_071415.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCC0
Protein Weight	61kDa
Shipping	Ice bag