


Anti-Mccc2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Hcc0

Item Number: CM0025965



Introduction

Rabbit polyclonal antibody against MCCC2 (MCCB) validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for leucine catabolism and mitochondrial carboxylase research.

[Learn More](#)

Parameter	Value
Product Name	MCCC2 Rabbit pAb
Remarks/Alias	MCCB; MCCCbeta; MCCC2
Species	Human
GeneID (Human)	64087
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 284-563 of human MCCC2 (NP_071415.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCC0
Protein Weight	61kDa
Shipping	Ice bag