

Anti-Mcee Rabbit Polyclonal Antibody For Wb, Elisa - Q96Pe7

Item Number: CM0018346



Introduction

Rabbit polyclonal antibody against human MCEE (methylmalonyl-CoA epimerase, mitochondrial), suitable for WB and ELISA. Cross-reacts with human and mouse. Research use for propionyl-CoA metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	MCEE Rabbit pAb
Remarks/Alias	MCE; MMCE; GLOD2; MCEE
Species	Human
Gene ID	84693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human MCEE (NP_115990.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96PE7
Protein Weight	19kDa
Shipping	Ice bag