

Anti-Mcl1 Rabbit Polyclonal Antibody For Wb - Q07820

Item Number: CM0023569



Introduction

Rabbit polyclonal antibody against human MCL1 (Mcl-1, Bcl2-L-3) for Western blot. Detects all isoforms; involved in apoptosis regulation. Peptide immunogen within amino acids 100-200.

[Learn More](#)

Parameter	Value
Product Name	MCL1 Rabbit pAb
Remarks/Alias	TM; EAT; MCL1L; MCL1S; Mcl-1; BCL2L3; MCL1-ES; bcl2-L-3; mcl1/EAT; MCL1
Species	Human
GeneID (human)	4170
GeneID	4170
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human MCL1 (NP_068779.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	Q07820
Protein Weight	37kDa
Shipping	Ice bag