

Anti-Mcm10 Rabbit Polyclonal Antibody For Wb And Elisa - Q7L590

Item Number: CM0024859

Introduction

Rabbit polyclonal antibody targeting human MCM10, a DNA replication initiation factor. Suitable for Western Blot and ELISA applications. Cross-reacts with human and mouse. Ideal for studying replication stress and genome stability.

[Learn More](#)

Parameter	Value
Product Name	MCM10 Rabbit pAb
Remarks/Alias	CNA43; DNA43; IMD80; PRO2249; MCM10
Species	Human
GeneID	55388
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 616-875 of human MCM10 (NP_877428.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7L590
Protein Weight	98kDa
Shipping	Ice bag