


Anti-Mcm2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49736

Item Number: CM0005776



Introduction

Rabbit monoclonal antibody against DNA replication licensing factor MCM2, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Suitable for studying cell cycle, DNA replication, and oncogenesis. Reliable IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	MCM2 Rabbit mAb
Remarks/Alias	BM28, CCNL1, CDCL1, cdc19, DFNA70, D3S3194, MITOTIN, MCM2
Species	Human
Gene ID	4171
Immunogen	Recombinant protein corresponding to amino acids 1-288 of human MCM2 (NP_004517.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49736
Protein Weight	102kDa
Shipping	Ice bag