


Anti-Mcm2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49736

Item Number: CM0002601



Introduction

Rabbit monoclonal antibody targeting human MCM2, validated for Western blot, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for DNA replication and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MCM2 Rabbit mAb
Remarks/Alias	BM28; CCNL1; CDCL1; cdc19; DFNA70; D3S3194; MITOTIN; MCM2
Species	Human
GeneID	4171
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MCM2 (P49736).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49736
Protein Weight	102kDa
Shipping	Ice bag