


Anti-Mcm2 Rabbit Polyclonal Antibody For Wb, Elisa - P49736

Item Number: CM0016309



Introduction

Rabbit polyclonal antibody targeting human MCM2, a DNA replication licensing factor. Validated for WB and ELISA, with cross-reactivity to human and mouse. Suitable for cell cycle and replication studies.

[Learn More](#)

Parameter	Value
Product Name	MCM2 Rabbit pAb
Remarks/Alias	BM28; CCNL1; CDCL1; cdc19; DFNA70; D3S3194; MITOTIN; MCM2
Species	Human
GeneID (Human)	4171
GeneID	4171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-700 of human MCM2 (NP_004517.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49736
Protein Weight	102kDa
Shipping	Ice bag