


Anti-Mcm3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P25205

Item Number: CM0003123



Introduction

High-quality rabbit monoclonal antibody against human MCM3 (P25205) for western blot, IHC, IF, IP, and ELISA. Validated in human, mouse, and rat samples. Ideal for DNA replication and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MCM3 Rabbit mAb
Remarks/Alias	HCC5; P1.h; RLFB; P1-MCM3; MCM3
Species	Human
GeneID (human)	4172
GeneID	4172
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25205
Protein Weight	91kDa
Shipping	Ice bag