


Anti-Mcm3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P25205

Item Number: CM0013573



Introduction

Rabbit polyclonal antibody targeting human MCM3. Validated for WB, IP, ELISA, and cross-reacts with human, mouse, rat. Suitable for DNA replication and cell cycle research. UniProt: P25205.

[Learn More](#)

Parameter	Value
Product Name	MCM3 Rabbit pAb
Remarks/Alias	HCC5; P1.h; RLFB; P1-MCM3; MCM3
Species	Human
GeneID (Human)	4172
GeneID	4172
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human MCM3 (NP_002379.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25205
Protein Weight	91kDa
Shipping	Ice bag