


Anti-Mcm5 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P33992

Item Number: CM0025689



Introduction

Rabbit polyclonal antibody against DNA replication licensing factor MCM5 (CDC46). Validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. For cell cycle and DNA replication studies.

[Learn More](#)

Parameter	Value
Product Name	MCM5 Rabbit pAb
Remarks/Alias	CDC46, MGORS8, P1-CDC46, MCM5
Species	Human
GeneID (Human)	4174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human MCM5 (NP_006730.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33992
Protein Weight	82kDa
Shipping	Ice bag