


Anti-Mcm5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P33992

Item Number: CM0014895



Introduction

Rabbit polyclonal antibody targeting human MCM5, validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat MCM5. Ideal for DNA replication and cell cycle studies.

[Learn More](#)

Parameter	Specification
Product Name	MCM5 Rabbit pAb
Remarks/Alias	CDC46; MGORS8; P1-CDC46; MCM5
Species	Human
GeneID (Human)	4174
GeneID	4174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human MCM5 (NP_006730.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33992
Protein Weight	82kDa
Shipping	Ice bag