


Anti-Mcm6 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Elisa - Q14566

Item Number: CM0002825



Introduction

KD validated anti-MCM6 rabbit monoclonal antibody for WB and ELISA; detects human and mouse MCM6 (Q14566), a key DNA replication licensing factor in the CMG helicase complex, ensuring once-per-cell-cycle replication.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MCM6 Rabbit mAb
Remarks/Alias	Mis5; P105MCM; MCG40308; [KO/KD Validated] MCM6
Species	Human
Gene ID (Human)	4175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 601-821 of human MCM6 (NP_005906.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14566
Protein Weight	93kDa
Shipping	Ice bag