


Anti-Mcm6 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q14566

Item Number: CM0020505



Introduction

Rabbit polyclonal antibody against human MCM6, validated for Western blot, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat. Ideal for DNA replication and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	MCM6 Rabbit pAb
Remarks/Alias	Mis5; P105MCM; MCG40308; MCM6
Species	Human
Gene ID (Human)	4175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 562-821 of human MCM6 (NP_005906.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14566
Protein Weight	93kDa
Shipping	Ice bag