


Anti-Mcmbp Polyclonal Antibody For Wb, Elisa - Q9Bte3

Item Number: CM0030999



Introduction

Anti-MCMBP Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human and mouse. Targets mini-chromosome maintenance complex-binding protein, a key regulator of DNA replication and sister chromatid cohesion. Ideal for cell cycle and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	MCMBP Rabbit pAb
Remarks/Alias	MCM-BP; C10orf119; MCMBP
Species	Human
Gene ID	79892
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human MCMBP (NP_079110.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BTE3
Protein Weight	73kDa
Shipping	Ice bag