


Anti-Mcoln2 (Trpml2) Rabbit Polyclonal Antibody For Wb And Elisa - Q8lzk6

Item Number: CM0018787



Introduction

Rabbit polyclonal antibody against human MCOLN2 (TRPML2), suitable for WB and ELISA. Cross-reacts with human and rat. Useful for studying lysosomal ion channel function, membrane trafficking, and immunity.

[Learn More](#)

Parameter	Value
Product Name	MCOLN2 Rabbit pAb
Remarks/Alias	TRPML2; TRP-ML2; MCOLN2
Species	Human
Gene ID (Human)	255231
Gene ID	255231
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-275 of human MCOLN2 (NP_694991.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8lzk6
Protein Weight	66kDa
Shipping	Ice bag