


Anti-Mct1/Slc16A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P53985

Item Number: CM0025699



Introduction

Rabbit polyclonal antibody against human MCT1/SLC16A1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat. Suitable for studying monocarboxylate transport and metabolism.

[Learn More](#)

Parameter	Value
Product Name	MCT1/SLC16A1 Rabbit pAb
Remarks/Alias	MCT; HHF7; MCT1; MCT1D; MCT1/SLC16A1
Species	Human
GeneID (Human)	6566
GeneID	6566
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 190-245 of human MCT1 (NP_003042.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53985
Protein Weight	54kDa
Shipping	Ice bag