

Anti-Mct4/Slc16A3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O15427

Item Number: CM0014117

Introduction

Rabbit polyclonal antibody targeting MCT4/SLC16A3, a monocarboxylate transporter. Validated for WB, IHC-P, ELISA in human, mouse, rat. Research use for lactate transport and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	MCT4/SLC16A3 Rabbit pAb
Remarks/Alias	MCT3; MCT4; MCT 3; MCT 4; MCT-3; MCT-4; MCT4/SLC16A3
Species	Human
GeneID (Human)	9123
GeneID	9123
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 410-465 of human MCT4/SLC16A3 (NP_004198.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15427
Protein Weight	49kDa
Shipping	Ice bag