

Anti-Mcu Monoclonal Antibody (Ko Validated) For Wb, Elisa - Q8Ne86

Item Number: CM0003439

Introduction

Rabbit monoclonal antibody against human MCU, validated for Western blot and ELISA. Recognizes human, mouse, and rat MCU (40 kDa). Suitable for studies on mitochondrial calcium uptake and cellular physiology.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MCU Rabbit mAb
Remarks/Alias	HsMCU; C10orf42; CCDC109A; CU
Species	Human
GeneID (Human)	90550
GeneID	90550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-170 of human MCU (NP_612366.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NE86
Protein Weight	40kDa
Shipping	Ice bag